

Smart-PSG



Bluetooth
Touchscreen
Connected

MADE IN FRANCE

**SLEEP
LABORATORY
EQUIPMENT**

ANALYSE

EFFICIENT

INTUITIVE

DIAGNOSE

INTERACTIVE

MOBILE

CONNECTED


RELIABLE

/ HOW IT WORKS



CID-LXr

- Placed **on the patient**,
- Measures **electrophysiological signals**,
- **Wireless communication** with TrackBoX,
- Used as a **polysomnograph** or a **polygraph**



TrackBoX

- **Real-time** display,
- **Programming** of recording,
- **Remote** supervision,
- **Touch** screen



TrackBoX Manager

- TrackBoX remote **visualisation**,
- Remote **management** of one or more Smart-PSG devices,
- **Supervision** of examinations



CIDELEC software

- Automatic signal **analysis**,
- Signal **interpretation**,
- **Data** archiving,
- **Customised reports**

These products are CE marked medical devices.

/ WE OFFER

CIDELEC provides:

- **Training and installation of products on site** by our team,
- **Innovative technology of our systems** to obtain precise and reliable measurements
- **Analysis and processing of information** with the creation of personalised reports,
- **After-sales service**, technical assistance and technical expertise



CIDELEC,
30 YEARS
YOUR PARTNER



Technical characteristics TrackBox

Dimensions: 290 x 215 x 155 mm- Weight: 3,2 kg

Channels	Bandwidth	Sampling frequency	Storage	Accuracy	Dots	Range
Pneumotachograph	0 - 10 Hz	4000 Hz	16 Hz	+/- 4 %	4096	+/- 1litre/s
Machine pressure	0 - 10 Hz	4000 Hz	256 Hz	+/- 25 Pa	4096	0 - 2 kPa
Analogue input	0 - 160 Hz	4000 Hz	256 Hz		24 bits = 16 777 215	+/- 1V
Luminosity		1 Hz	1 Hz		255	0 - 100 Lux



Technical characteristics CID-LXr

Dimensions: 32 x 82 x 114 mm - Weight: 152 g- Battery: Li-Po 1700 mAh - 3,7 V

Channels	Bandwidth	Sampling frequency	Storage	Accuracy	Dots	Range	Others
Breathing sound	200 - 2000 Hz	4000 Hz	Sound intensity at16 Hz		256		Sensitivity 20 - 80 dB Adaptive threshold
Snoring sound	20 - 200 Hz	4000 Hz	Sound intensity at 16 Hz		256		Sensitivity 60 - 120 dB Adaptive threshold
Supra-sternal pressure	0,02 - 20 Hz	4000 Hz	8 Hz		4096	+/- 100 Pa	
Position		1 Hz	1 Hz				5 positions
Actimetry		1000 Hz	8 Hz				
Nasal flow	0 - 10 Hz	4000 Hz	256 Hz		65536	+/- 300 Pa	
SpO ₂ ⁽¹⁾			1 Hz	+/- 3 % (from 70 to 100 %) ⁽²⁾	100	0 - 100 %	Averaged over 4 pulse cycles
Pulse frequency ⁽¹⁾			1 Hz	+/- 5 BPM ⁽²⁾	256	40 - 240 BPM ⁽²⁾	
Photoplethysmogram ⁽¹⁾			64 Hz				
Inductive belts	0,1 - 10 Hz	32 Hz	8 Hz		65536		
6 EEG channels	45 Hz	500 Hz	250Hz		65536	1627 µV	Built-in 50 Hz noise tester
2 EOG channels	45 Hz	500 Hz	250Hz		65536	1627 µV	Built-in 50 Hz noise tester
1 ECG channel	45 Hz	500 Hz	250Hz		65536	6510 µV	Built-in 50 Hz noise tester
4 EMG channels	10 - 100 Hz	4000 Hz	64Hz		256		

(1) Nonin manufacturer
(2) Under worst-case conditions



Take advantage
of CIDELEC
technology with
PneaVoX
sensor accuracy



Remote supervision
Bluetooth connectivity
Signal analysis

OUR OBJECTIVE : ENABLING YOU TO USE ACCURATE,
COMPLETE AND RELIABLE RECORDINGS.

/ PNEAVOX

PneaVoX technology is unique.

One sensor records 3 physiological parameters :

- Buccal and nasal **breathing**,
- **Respiratory effort** via suprasternal pressure to differentiate between obstructive, central and combined apneas,
- **Snoring** (energy, intensity).

Finally, the PneaVoX sound sensor **analyses upper airway resistance** by measuring the sound intensity.

«The **PneaVoX sound sensor**, to improve differentiation between sleep disorders via the analysis of tracheal sounds.»



/ SCIENTIFIC BIBLIOGRAPHY

M. Glos, A. Sabil, K.S. Jelavic, C; Schöbel, I. Fietze, T. Penzel.
Characterization of respiratory events in obstructive sleep apnea using suprasternal pressure monitoring.
J Clin Sleep Med. 2018; 14(3): 359-369.

Penzel T, Sabil A. ***The use of tracheal sounds for the diagnosis of sleep apnoea.***
Breathe 2017; 13: e37-e45.

A. Amadeo, M. Fernandez-Bolanos, J.O. Arroyo, S. Khirani, G. Baffet, B. Fauroux. ***Validation of a Suprasternal Pressure Sensor for Sleep Apnea Classification in Children,***
Journal of Clinical Sleep Medicine, Vol. 12, No. 12, 2016.

[...]



20 rue des Métiers

49130 Sainte Gemmes sur Loire - FRANCE

Tel : +33 (0)2 41 66 20 88 - Fax : +33 (0)2 41 79 07 76

Email : sales@cidelec.net

Website : www.cidelec.net

The content of this brochure (text, illustrations, schemes, images, etc...) is the exclusive property of CIDELEC.
The CIDELEC logo and the CIDELEC, CID, CID102, TrackBoX and PneaVoX brands are registered by CIDELEC.
CIDELEC reserves the right to modify this document at anytime.

The SMART-PSG is a medical device classified IIa range, manufactured by the CIDELEC company - CE N°0459
The SMART-PSG is a system for recording physiological signals, for the diagnosis of sleep disorders.

Read the instructions carefully before use.

Document modified in 07/2023